

PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public.class-woocommerce-print-products-public.php on line 2806 bMC275 Clone ID

HSA; PRO0883; PRO0903; PRO1341 Synonyme

Host Species

Rabbit Anti-ALB antibody(DMC275); IgG1 Chimeric mAb Description

Delivery In Stock

Uniprot ID P02768

Rabbit/Human Fc chimeric IgG1 IgG type

Clonality Monoclonal Reactivity Human Application Flow Cyt Recommend Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Score at 20°C to 9.0°C for 12 months in yoshilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid Score at 20°C to 9.0°C for 12 months in yoshilized proteins are shipped at ambient temperature.

The gene free for the start of the start Storage & Shipping

Background

Usage Conjugate Research use only Unconjugated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.




Figure 1. Flow cytometry analysis with Anti-ALB (DMC275) on HEK293 cells transfected with human ALB (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

Email: info@dimabio.com Website: www.dimabio.com

