

PRODUCT INFORMATION

| | |
|-----------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Common Name | Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\publicclass\woocommerce-print-products-public.php on line 2806 LY:3372993;LY:3372993;LY3372993;N3pG 4;N3pG4;N3PG-IV |
| Conjugate | Unconjugated |
| Synonyms | AAA;AD1;PN2;ABPP;APPI;CVAP;ABETA;PN-II;preA4;CTFgamma;alpha-sAPP |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | N/A |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | APP |
| Uniprot ID | P05067 |
| Description | Anti-APP(remternetug biosimilar)(IgG1) mAb |
| Delivery | 4-6 weeks |
| Yefei Storage | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

DIMABIO CONFIDENTIAL

