

PRODUCT INFORMATION

| | |
|---|--|
| Common Name | Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 CH 7092230;CH7092230;CH-7092230;RG 6107;RG-6107;RG6107;RO 7092230;RO-7112689;RO7112689;SKY 59;SKY-59;SKY59 |
| Conjugate | Unconjugated |
| Synonyms | C5D;C5a;C5b;ECLZB;CPAMD4 |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | N/A |
| IgG type | IgG1 |
| Reactivity | Human |
| Target | C5 |
| Uniprot ID | P01031 |
| Description | Anti-C5(crovalimab biosimilar[IgG1]) mAb |
| Delivery | 4-6 weeks |
| Yefei Storage | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

DIMABIO CONFIDENTIAL

