

## PRODUCT INFORMATION

<b>Common Name</b>	<b>Warning:</b> Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 MOR208, XmAb5574, XmAb@5574, MOR-00208, MOR-208, XENP-5574, XMP-5574, XmAb-5574, Xmab-CD19, Xmab5574, tafasitamab-cxix
<b>Conjugate</b>	Unconjugated
<b>Synonyms</b>	CD19,B4,CVID3,MGC12802
<b>Applications</b>	ELISA
<b>Recommended Dilutions</b>	ELISA 1:5000-10000
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Host Species</b>	Humanized
<b>IgG type</b>	Human IgG1-G2 - kappa
<b>Reactivity</b>	Human
<b>Target</b>	CD19
<b>Uniprot ID</b>	P15391
<b>Description</b>	Anti-CD19(tafasitamab biosimilar) mAb
<b>Delivery</b>	In Stock
<b>Yefei Storage</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only

## Anti-CD19(tafasitamab biosimilar) mAb ELISA

0.2 µg of Human CD19 (M75I,L82V,F83L), hFc tagged protein per well




Figure 1. ELISA plate pre-coated by 2 µg/mL (100 µL/well) Human CD19(M75I,L82V,F83L), hFc Tag (PME101557) can bind Anti-CD19(tafasitamab biosimilar) mAb (BME100683) in a linear range of 3.2-80 ng/mL.

