

PRODUCT INFORMATION

Warning: Undefined variable \$hacAttributeValueDescription in CAywwwroot\mirror.dlmablo.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 311 Clone ID

CD200R1

CD200R;HCRTR2;MOX2R;OX2R Synonyme

Host Species Rabbit

Anti-CD200R1 antibody(3H1), IgG1 Chimeric mAb Description

Delivery In Stock Uniprot ID Q8TD46

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Applicatio Human Flow Cyt Recommend Dilutions Flow Cyt 1/100

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

specific instructions or reconstitution. Some at 20°C to 80°C for 12°C to 80°C for 12°C to 80°C for 12°C months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at a mishent temperature. This gene encodes a receptor for the OX.2 membrane glycoprotein. Both the receptor an substrate are cell surface glycoproteins containing two This gene encodes a receptor for the OX.2 membrane glycoproteins. Both the receptor and substrate are cell surface glycoproteins containing two This gene encodes are receptor and the substrate in the substrate interaction may function as a myeloid downregulatory signal. Mouse studies of a related gene suggest that this interaction may control myeloid function in a tissue-specific manner. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq. Jul 2008] Background

Usage Conjugate Research use only Unconjugated




Figure 1. Flow cytometry analysis with $10\mu g/mL$ Anti-CD200R1 (3H1) mAb on HEK293 cells transfected with human CD200R1 (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

Email: info@dimabio.com Website: www.dimabio.com

