

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 COR-004;COR004;COR 004;COR-4;COR4;COR 4;CPI-006;CPI006;CPI 006;CPI-6;CPI6;CPI 6;CPX-006;CPX006;CPX 006;CPX-6;CPX6;CPX 6
Conjugate	Unconjugated
Synonyms	NT;eN;NT5;NTE:eNT;CD73;E5NT;CALJA
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG1
Reactivity	Human
Target	CD73
Uniprot ID	P21589
Description	Anti-CD73(mupadolimab biosimilar)(IgG1) mAb
Delivery	4-6 weeks
Yefei Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

DIMABIO CONFIDENTIAL

