

## **PRODUCT INFORMATION**

Warning . Undefined variable \$hasAttributeValueDescription in Calwawroottmirror.dimabio.com/wp-content\plugins/woocommerce-print-products-public.php on line 2806 Chongqing Precision

Conjugate

Unconjugated
CEACAM-5,CD66e,CEA,Meconium antigen 100 Synonyms

Applications Flow Cyt Recommend Dilutions Flow Cyt 1:100

Target

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. Formulation & Reconstitution

IgG type Human IgG1 - kappa Reactivity

Uniprot ID Anti-CEACAM5 (Chongging Precision biosimilar) mAb Description

CEACAM5

Delivery

In Stock
Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals

Usage




Figure 1. Flow cytometry analysis with  $1\mu$ g/mL Anti-CEACAM5 (Chongqing Precision biosimilar) mAb (BME100182) on HEK293 cells transfected with Human CEACAM5 protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

Email: info@dimabio.com Website: www.dimabio.com

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

