

**PRODUCT INFORMATION**

<b>Common Name</b>	AMG-757 without CD3
<b>Synonyms</b>	SCDO1
<b>Applications</b>	Flow Cyt
<b>Recommended Dilutions</b>	Flow Cyt 1:100
<b>Formulation &amp; Reconstitution</b>	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
<b>Host Species</b>	Humanized
<b>IgG type</b>	IgG1-kappa
<b>Reactivity</b>	Human
<b>Target</b>	DLL3
<b>Uniprot ID</b>	Q9NYJ7
<b>Description</b>	Anti-DLL3(tarlatamab without CD3 biosimilar) mAb
<b>Delivery</b>	In Stock
<b>Storage &amp; Shipping</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
<b>Background</b>	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated
<b>DIMA Disclaimer</b>	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.






Figure 1. Flow cytometry analysis with 15µg/mL Anti-DLL3(tarlatamab without CD3 biosimilar) mAb(BME100257) on 293T-DLL3 cells (Blue histogram) or isotype control mAb on 293T-DLL3 cells (Red histogram).

