

PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public:class-woocommerce-print-products-public-php on line 2806 bM95 Clone ID

DNAM1

DNAM1; CD226; PTA1 Synonyme

Host Species Rabbit

Anti-DNAM-1 antibody(DM95); Rabbit mAb Description

Delivery In Stock Uniprot ID Q15762

IgG type Rabbit IgG Clonality Monoclonal Reactivity Human Application ELISA; Flow Cyt

Recommend Dilutions ELISA 1:5000-10000; Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Scre at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperated freezing and thawing). Lyophilized proteins are shipped at ambient temperated and subset of T cells. It is a member of the lg-superfamily containing 2 ig-like domains of the V-set. The protein mediates cellular adhesion of platelets and megacycyctic cells to vascular endothelial cells. The protein also plays a role in megakaryocycic reliad maturation. Atternative spiking results in multiple transcript variants. Storage & Shipping

Usage Research use only

Necessity used only

Unconjugated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer




Figure 1. ELISA plate pre-coated by 2 μg/ml (100 μl/well) Human DNAM-1 protein, mFc-His tagged protein ([getskuurl sku="PME100050"]) can bind Rabbít anti-DNAM-1 monoclonal antibody (clone: DM95) in a linear range of 0.64-80 ng/ml.

Email: info@dimabio.com Website: www.dimabio.com

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)









Figure 2. HEK293 cell line transfected with irrelevant protein **(A)** and human DNAM-1 **(B)** were surface stained with Rabbit anti-DNAM-1 monoclonal antibody 1µg/ml (**clone: DM95**) followed by PE-conjugated anti-rabbit IgG secondary antibody.



Email: info@dimabio.com Website: www.dimabio.com

