

PRODUCT INFORMATION

Warning: Undefined variable \$hasAttributeValueDescription in C\\u00e4wwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products-public.php on line 2806
Z84 Clone ID

Host Species Rabbit

Anti-Eribulin antibody(2E4); Rabbit mAb Description

Delivery In Stock Uniprot ID N/A lgG type Rabbit IgG Clonality Monoclonal Reactivity ELISA Applicati

Recommend Dilutions ELISA 1/5000-10000

Purification Purified from cell culture supernatant by affinity chromatography

Store at 20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

N/A Formulation & Reconstitution

Storage & Shipping

Research use only Usage Conjugate

Unconjugated
All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaime

ELISA assay to evaluate Eribulin-Antibody 0.2µg Human IgG-Eribulin per well




Figure 1. Elisa plates were pre-coated with IgG-Eribulin (0.2 μ g/per well). Serial diluted anti-Eribulin monoclonal antibody (DME101047) solutions were added, washed, and incubated with secondary antibody before Elisa reading. From above data, the EC50 for anti-Eribulin monoclonal antibody binding with IgG-Eribulin is 1.658ng/ml.



Email: info@dimabio.com Website: www.dimabio.com

