

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttribute\$ in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\publicclass-woocommerce-print-products-public.php on line 2806 NVS250519
Conjugate	Unconjugated
Synonyms	FA11; FXI; PTA
Applications	ELISA
Recommended Dilutions	ELISA 1:5000-10000
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Homo sapiens
IgG type	Human IgG1(K214R,D265A,P329A) - lambda2
Reactivity	Human
Target	F11
Uniprot ID	P03951
Description	Anti-F11(abelacimab biosimilar) mAb
Delivery	In Stock
Yefei Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

Anti-F11(abelacimab biosimilar) mAb ELISA

0.2 µg of Human F11, FLAG tagged protein per well




Figure 1. ELISA plate pre-coated by 2 µg/mL (100 µL/well) Human F11 Protein, Flag Tag(PME101748) can bind Anti-F11(abelacimab biosimilar) mAb (BME100255) in a linear range of 0.13-16 ng/mL.

