

PRODUCT INFORMATION

Common Name	M2
Synonyms	DDK□Flag□
Applications	ELISA, Flow Cyt
Recommended Dilutions	ELISA 1:5000-10000, Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Mouse
IgG type	IgG1-kappa
Reactivity	N/A
Target	Flag
Uniprot ID	N/A
Description	Anti-Flag (BM301) Mouse mAb
Delivery	In Stock
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).Lyophilized antibodies are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	Unconjugated
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.



Anti-Flag (BM301) Mouse mAb ELISA

0.2 μg of Human LRP6(20-1370), Flag tagged protein per well




Figure 1. ELISA plate pre-coated by 2 $\mu\text{g}/\text{mL}$ (100 $\mu\text{L}/\text{well}$) Human LRP6(20-1370) Protein, Flag Tag (PME101719) can bind Anti-Flag (BM301) Mouse mAb (BME100093) in a linear range of 0.64-80 ng/mL.

