

PRODUCT INFORMATION

Common Name	AMG-477;AMG477;AMB 477;REMD-477;REMD477;REMD 477
Synonyms	GR;GCR;GRL;GCCR;GCRST
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG2
Reactivity	Human
Target	GCCR
Uniprot ID	P04150
Description	Anti-GCCR(volagidemab biosimilar) IgG2) mAb
Delivery	3-4 weeks
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	Unconjugated
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

DIMA BIOTECH CONFIDENTIAL

