

PRODUCT INFORMATION

Warning: Undefined variable \$hacAttributeValueDescription in CAywwwroot\mirror.dlmablo.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 Clone ID

GLP1R

GLP-1;GLP-1-R;GLP-1R

Host Species Rabbit

Anti-GLP1R antibody(1A4), Rabbit mAb Description In Stock

Delivery Uniprot ID P43220 IgG type Rabbit IgG Clonality Monoclonal Reactivity Human Flow Cyt Application Recommend Dilutions Flow Cyt 1/100

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

Background

specific instructions of reconstitution.

Store at 2-0°C to 8-0°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

This gene encodes a 7-transmembrane protein that functions as a receptor for glucagon-like peptide 1 (GLP-1) hormone, which stimulates glucose-induced insulin secretion. This receptor, which functions at the cell surface, becomes internalized in response to GLP-1 and GLP-1 analogs, and it plays an important role in the signaling cascades leading to insulin secretion. It also displays neuroprotective effects in alm models. Polymorphisms in this gene are associated with diabetes. The protein is an important drug target for the treatment of type 2 diabetes and stroke. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Apr 2016]

Research use only

Conjugate Unconjugated




Figure 1. Flow cytometry analysis with 2µg/mL Anti-GLP1R (1A4) mAb on HEK293 cells transfected with human GLP1R (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

Email: info@dimabio.com Website: www.dimabio.com



