

## **PRODUCT INFORMATION**

Warning: Undefined variable shasAttributeValueDescription in C\u00e4wwroot\u00e4mirror.dimabio.com\u00fcp-content\plugins\u00e4woocommerce-print-products\u00e4public\u00e4class-woocommerce-print-products-public.php on line 2806 612 Clone ID

Human IgG Synonyme Host Species Rabbit

Anti-Human IgG antibody(6E12), Rabbit mAb Description Delivery In Stock

Uniprot ID N/A IgG type Clonality Rabbit IgG Monoclonal Reactivity
Applications Human ELISA

Recommended Dilutions ELISA 1:5000-10000

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Specific in BioCoccusts of reconstructions. Store at 2-20°C to 80°C (are 12°C to 80° Storage & Shipping

Usage Research use only Conjugate

Unconjugated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

SKU	Clone ID	Species Reactivity			Cross Reactivity		
		Mouse IgG	Rat IgG	Goat IgG	Human-IgG1	Human-IgG2	Human-IgG4
DME101028	1G9	92	( ×	<u>201</u>	+ + +	+	+ + +
DME101029	2D4	-	+	ā	+ +	+ +	-
DME101030	2D9	4	+ + +	8	+ + +	+ +	-
DME101031	2H2	22	=	12	-	=	22
DME101032	3D11	+	+	+	+ + +	+ +	+ + +
DME101033	4A12	<b>+</b>	+	+	-	=	) <del></del>
DME101034	5B4	+	+	+	+ + +	+ +	+++
DME101037	6C5	92	2	355	+ + +	+ +	+ + +
DME101035	6E12		2	la la	+ + +	+ +	+ + +
DME101036	1E12	+	+	+	+ + +	+	+ + +
DME101018	1D9	12	Ψ,	≅	+ + +	+	+ + +

Figure 1. ELISA examination of rabbit anti-human IgG mAbs binding to immunogolobulins of different species and isotypes.

Email: info@dimabio.com Website: www.dimabio.com

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

