Cat. No. BME100535



PRODUCT INFORMATION

Warning: Undefined variable \$hasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products-public.php on line 2806 AU-007;AU-007;BD-8;BD8;BD 8;BD68;BD G8;BD G8 Common Name

Conjugate

Unconjugated p55;CD25;IL2R;IMD41;TCGFR;IDDM10;IL2-RA;IL-2-RA Synonyms

Applications Recommended Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species lgG1 IgG type Reactivity Target IL2RA Uniprot ID

Anti-IL2RA(imneskibart biosimilar||IqG1) mAb Description

Delivery

3.4 weeks
Store at .20°C to .80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at .80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

> Email: info@dimabio.com Website: www.dimabio.com

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Usage

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

