

## **PRODUCT INFORMATION**

Warning. Undefined variable \$hasAttributeValueDescription in Calwawroottmirror.dimablo.com/wp-content\plugins/woocommerce-print-products-public.php on line 2806
8A7-1905284.

Conjugate Unconjugated C1orf32; Angulin-3 Synonyms Applications Flow Cyt

Recommend Dilutions Flow Cyt 1:100 Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Human IgG2(S378A) - Kappa IgG type

Reactivity Target ILDR2 Uniprot ID

Description Anti-ILDR2(bapotulimab biosimilar) mAb

Delivery

In Stock
Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals

Usage




Figure 1. Flow cytometry analysis with  $1\mu g/mL$  Anti-ILDR2(bapotulimab biosimilar) mAb (BME100243) on HEK293 cells transfected with Human ILDR2 protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).



Email: info@dimabio.com Website: www.dimabio.com

