

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttribute\$Value\$Description in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\publicclass-woocommerce-print-products-public.php on line 2806 BAY-1905254
Conjugate	Unconjugated
Synonyms	C1orf32; Angulin-3
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Homo sapiens
IgG type	Human IgG2(S378A) - Kappa
Reactivity	Human
Target	ILDR2
Uniprot ID	Q71H61
Description	Anti-ILDR2(bapotulimab biosimilar) mAb
Delivery	In Stock
Yefei Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only




Figure 1. Flow cytometry analysis with 1 μ g/mL Anti-ILDR2(bapotulimab biosimilar) mAb (BME100243) on HEK293 cells transfected with Human ILDR2 protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

