

PRODUCT INFORMATION

Synonyms

Unconjugated Conjugate Recommended Dilutions Flow Cyt 1:100

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species Humanized IgG type Reactivity Human IgG2 - Kappa Human Target Uniprot ID ITGAV P06756

Description Anti-ITGAV(abituzumab biosimilar) mAb Delivery In Stock

Store at 20°C to .80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Research use only Storage & Shipping

Background

Usage




Figure 1. ITGAV protein is highly expressed on the surface of A549 cell membrane. Flow cytometry analysis with 1 μ g/mL Anti-ITGAV(abituzumab biosimilar) mAb (BME100146) (Blue histogram) or isotype control mAb (Red histogram) on A549





