

PRODUCT INFORMATION

Common Name	SGN B6A□Sigvotatug
Synonyms	Integrin beta-6
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Humanized
IgG type	IgG1(E239D,M241L)- kappa
Reactivity	Human
Target	ITGB6
Uniprot ID	P18564
Description	Anti-ITGB6(sigvotatug biosimilar) mAb
Delivery	In Stock
Storage & Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only
Conjugate	Unconjugated
DIMA Disclaimer	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

CONFIDENTIAL






Figure 1. Flow cytometry analysis with 1µg/mL Anti-ITGB6(sigvotatug biosimilar) mAb(BME100248) on Expi293 cells transfected with Human ITGB6&ITGAV protein (Blue histogram) or Expi293 transfected with Human ITGAV protein (Red histogram).

DIMABIO CONFIDENTIAL

