

**PRODUCT INFORMATION**

<b>Clone ID</b>	9E10
<b>Target</b>	Myc
<b>Synonyms</b>	9E10 antibody; 9E10 anti myc; 9E10 myc; anti-Myc; anti-myc antibodies; Anti-myc Antibody; anti c-myc; anti cmc; anti myc; c-myc antibody; cmc antibody; myc antibody; tag antibody
<b>Host Species</b>	Mouse
<b>Description</b>	Anti-Myc Monoclonal Antibody
<b>Delivery</b>	In Stock
<b>Uniprot ID</b>	N/A
<b>IgG type</b>	Mouse IgG
<b>Clonality</b>	Monoclonal
<b>Reactivity</b>	Myc tag
<b>Applications</b>	FC, IF, IP, WB
<b>Recommended Dilutions</b>	WB: 1:500-1:2000, FC: 1:100, IF 1:100, IP 1:100-200
<b>Purification</b>	N/A
<b>Formulation &amp; Reconstitution</b>	0.64 mg/ml in PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide
<b>Storage &amp; Shipping</b>	Stable for 1 year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to opening the cap. Aliquot to avoid repeated freezing and thawing. Aliquot will be stable at 4°C for 3 months.
<b>Background</b>	N/A
<b>Usage</b>	Research use only
<b>Conjugate</b>	Unconjugated
<b>DIMA Disclaimer</b>	All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

