

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttribute\$ in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\publicclass-woocommerce-print-products-public.php on line 2806 NP-137
Conjugate	Unconjugated
Synonyms	NTN1; NTN1L
Applications	ELISA
Recommended Dilutions	ELISA 1:5000-10000
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	Netrin-1
Uniprot ID	Q95631
Description	Anti-Netrin-1(NP-137 biosimilar) mAb
Delivery	In Stock
Yefei Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

Anti-Netrin-1(NP-137 biosimilar) mAb ELISA

0.2 µg of Human Netrin-1, Flag tagged protein per well




Figure 1. ELISA plate pre-coated by 2 µg/mL (100 µL/well) Human Netrin-1 Protein, Flag Tag (PME101766) can bind Anti-Netrin-1(NP-137 biosimilar) mAb (BME100199) in a linear range of 0.64–16 ng/mL.

