

PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\u00e4wwroot\u00fcmirror.dimabio.com\u00fcwp-content\plugins\u00fcwoocommerce-print-products-public.php on line 2806
AGEN:2034-AGEN2034;Remba,7-00;Remba70;0-00; Common Name

Conjugate

Synonyms

Applications Recommended Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species IgG4 IgG type Reactivity Target PD-1 Uniprot ID

Anti-PD-1(balstilimab biosimilar∏lgG4) mAb Description

Delivery

3.4 weeks
Store at .20°C to .80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at .80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

> Email: info@dimabio.com Website: www.dimabio.com

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Usage

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

