

PRODUCT INFORMATION

Warning: Undefined variable \$hasAttributeValueDescription in C:\www.root\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products-public.php on line 2806
HXX-10HX10P.HX10 Common Name

Conjugate

PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1 Synonyms

Applications Recommended Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution Host Species

IgG4 IgG type Reactivity Target PD-1 Uniprot ID

Anti-PD-1(serplulimab biosimilar∏lqG4) mAb Description

Delivery

A-4 weeks
Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

> Email: info@dimabio.com Website: www.dimabio.com

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Usage

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

