

PRODUCT INFORMATION

Common Name

Conjugate

Synonyms Applications PD1;PD-1;CD279;SLEB2;hPD-1;hPD-I;hSLE1

Recommended Dilutions N/A

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Host Species N/A IgG4 Reactivity Human Target PD-1 Uniprot ID Q15116

Anti-PD-1(zimberelimab biosimilar[]lgG4) mAb

Delivery 3-4 weeks

See at 20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature than the store at -80°C (Avoid Research grade biosimilar. Not for use in therapeut or diagnostic procedures from humans or animals. Storage & Shipping

> Email: info@dimabio.com Website: www.dimabio.com

Background Usage Research use only



