

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 Cejemly;CS-1001;CS1001;CS 1001;Sugemalimab;WBP-3155;WBP3155;WBP 3155;Zegema;Zegemet;Zegemex;ZEJEME1 |
| Conjugate | Unconjugated |
| Synonyms | B7-H;B7H1;PDL1;PD-L1;hPD-L1;PDCD1L1;PDCD1LG1 |
| Applications | N/A |
| Recommended Dilutions | N/A |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | N/A |
| IgG type | IgG4 |
| Reactivity | Human |
| Target | PD-L1 |
| Uniprot ID | Q9NZQ7 |
| Description | Anti-PD-L1(sugemalimab biosimilar)(IgG4) mAb |
| Delivery | 4-6 weeks |
| Yefei Storage | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |

DIMABIO CONFIDENTIAL

