

PRODUCT INFORMATION

Warning: Undefined variable \$hasAttributeValueDescription in C\u00e4wwroot\u00fanirror.dimabio.com\u00e4wp-content\plugins\u00e4woocommerce-printproducts\u00e4public\u00e4class-woocommerce-print-products-public.php on line 2806
487-Clone ID

THEMIS

THEMIS, C6orf190, C6orf207 Synonyme

Host Species Rabbit

Anti-pTHEMIS antibody(4F2), IgG1 Chimeric mAb Description

Delivery In Stock Uniprot ID Q8N1K5

Rabbit/Human Fc chimeric IgG1 lgG type

Clonality Monoclonal Reactivity Human ELISA Application

Recommend Dilutions ELISA 1:5000-10000

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

This gene encodes a protein that plays a regulatory role in both positive and negative T-cell selection during late thymocyte development. The protein functions through T-cell antigen receptor signaling, and is necessary for proper lineage commitment and maturation of T-cells. Alternative splicing results in multiplie transcript variants. Storage & Shipping

Usage Research use only

Inconjugated
All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer




Figure 1. ELISA examination of DMC101054 specificity. The antibody at various concentrations was incubated in plates coated with either THEMIS or phosphorylated THEMIS (pTHEMIS), and a sigmoidal curve was developed only when using pTHEMIS, but not THEMIS, as capture antigen, indicating a specific binding of DMC101054 to pTHEMIS.

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

Email: info@dimabio.com Website: www.dimabio.com

