

PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\u00e4wroot\u00fcmirror.dimablo.com\u00e4wp-content\plugins\u00e4woocommerce-print-products\u00e4public\u00e4class-woocommerce-print-products-public.php on line 2806
bM25 Clone ID

S protein SARS-CoV-2 Spike S2 Synonyme

Host Species Rabbit

Anti-SARS-CoV-2 Spike antibody(DM25); Rabbit mAb Description Delivery In Stock

Uniprot ID PODTC2 lgG type Rabbit IgG Clonality Monoclona Reactivity SARS-CoV-2 ELISA Application

Recommend Dilutions ELISA 1:5000-10000

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis fo specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

Score at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at another temperature. As a strength of the stren Background

Usage Conjugate Research use only Unconjugated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement.

Anti-SARS-CoV-2 Spike S2 Antibody Elisa

0.2 µg of SARS-CoV-2 Spike S2 protein per well




Figure 1. Elisa plate pre-coated by 2 μg/ml(100μl/well) SARS-CoV-2 Spike S2 protein can bind Rabbit Anti-SARS-CoV-2 Spike S2 monoclonal antibody (clone:DM25) in a linear range of 1.9-250 ng/ml.

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

Email: info@dimabio.com Website: www.dimabio.com

