Cat. No. BME100427



PRODUCT INFORMATION

Warning: Undefined variable \$hasAttributeValueDescription in C:\www.root\mirror.products\public\class-woocommerce-print-products-public.php on line 2806 EOS 884448;EOS-88448;EOS-849448;EOS-849448;EOS-849448;EOS-849468;EOS-8494895;EOS-8496;EOS-849 Common Name

Unconjugated VSIG9;VSTM3;WUCAM Conjugate Synonyms

Applications N/A Recommended Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species

lgG1 IgG type Reactivity Human Target TIGIT Uniprot ID

Anti-TIGIT(belrestotug biosimilar∏lgG1) mAb Description

Delivery

3.4 weeks
Store at .20°C to .80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at .80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

> Email: info@dimabio.com Website: www.dimabio.com

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Usage

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

