

PRODUCT INFORMATION

Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 DPR-305

Unconjugated CD282;TIL4 Conjugate Synonyms Applications Flow Cyt Recommende Dilutions Flow Cyt 1:100

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species

Human IgG4 - Kappa IgG type Reactivity Target TLR2 Uniprot ID 060603

Anti-TLR2(tomaralimab biosimilar) mAb Description

Delivery

In Stock
Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shi

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals

Usage




Figure 1. Flow cytometry analysis with $15\mu g/mL$ Anti-TLR2(tomaralimab biosimilar) mAb (BME100110) on HEK293 cells transfected with Human TLR2 protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

Email: info@dimabio.com Website: www.dimabio.com

