

PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C\\u00edwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public.php on line 2806
b5 1062a\u00edb1062a\u00edb5-1062a\u00edb1062a\u00edb5-1062a\u00edb1062a\u00edb5-1062a\u00edb1062a\u00edb5-1062a\u00edb1062a\u00edb5-1062a\u00edb1062a\u00edb5-1062a\u00edb1062a\u00edb5-1062a\u00edb1062a\u00edb5-1062a\u00edb10 Common Name

Conjugate

Unconjugated EGP1;GP50;M1S1;EGP-1;TROP2;GA7331;GA733-1;TACSTD2 Synonyms

Applications Recommended Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species lgG1 IgG type Reactivity Target Trop2 Uniprot ID

Anti-Trop2(datopotamab biosimilar∏lqG1) mAb Description

Delivery

3.4 weeks
Store at .20°C to .80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at .80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

> Email: info@dimabio.com Website: www.dimabio.com

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Usage

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

