

PRODUCT INFORMATION

\$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-commerce-print-products-public.php on line 2806

Conjugate

Unconjugated
TNFRSF12A;FGF-inducible 14;FN14;TweakR;CD266 Synonyms

Applications ELISA; Flow Cyt

Recommend Dilutions ELISA 1:5000-10000; Flow Cyt 1:100

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

IgG type Human IgG1 - kappa Reactivity Target TweakR Uniprot ID

Description Anti-TweakR (enavatuzumab biosimilar) mAb

Delivery

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals

Usage

Anti-TweakR (enavatuzumab biosimilar) mAb ELISA




Figure 1. ELISA plate pre-coated by 2 μ g/ml (100 μ l/well) Human TweakR, mFc tagged protein ([getskuurl sku="PME100463"]) can bind Anti-TweakR Neutralizing antibody in a linear range of 0.64-80 ng/ml.



Email: info@dimabio.com Website: www.dimabio.com






Figure 2. Flow cytometry analysis with 1 μ g/mL Anti-TweakR (enavatuzumab) mAb (BME100031) on HEK293 cells transfected with Human TweakR protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

Email: info@dimabio.com Website: www.dimabio.com

