

PRODUCT INFORMATION

shasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-ommerce-print-products-public.php on line 2806 Clone ID

HEL-117; HEL-S-53; NDGOA; PARK5; PGP 9.5; PGP9.5; PGP95; SPG79; Uch-L1

Host Species

Anti-UCHL1 antibody(2A4), Rabbit mAb Description

Delivery In Stock Uniprot ID P09936 lgG type Rabbit IgG Clonality Monoclonal Reactivity ELISA; WB Applicati

Recommend Dilutions ELISA 1:5000-10000; WB 1:1000

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

Scre at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. The protein encoded by this gene belongs to the peptidase CLT family. This enzyme is a thiol protease that hydrolyzes a peptide bond at the C-terminal glycine of ubiquitin. This gene is specifically expressed in the neurons and in cells of the diffuse neuroendocrine system. Mutations in this gene may be associated with Parkinson disease, [provided by ReFegs. Spe 2009]

Usage Research use only

Unconjugated
All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any proteir
engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

Anti-UCHL1 (2A4) mAb ELISA 0.1 µg of Human UCHL1, hFc tagged protein per well




Figure 1. ELISA plate pre-coated by 1 μ g/ml (100 μ l/well) Human UCHL1 protein, hFc tagged protein (PME100629) can bind Rabbit anti-UCHL1 monoclonal antibody(clone: 2A4) in a linear range of 1-20 ng/ml.

Email: info@dimabio.com Website: www.dimabio.com









Figure 2.Anti-UCHL1 antibody (SKU# DME100843) at 1/1000 dilution

Lane: RAJI(human Burkitt's lymphoma B lymphocyte), whole cell lysate

Email: info@dimabio.com Website: www.dimabio.com

Secondary: Goat Anti-Rabbit IgG H&L (HRP) at 1/5000 dilution

band size: 24 kDa

