

PRODUCT INFORMATION

Target	AXL, CD63, EPCAM, LGALS1
Description	The embedded and sliced tumor sections (4µm section) from Balb/C nu mice implanted with PC3 cells
Host Cells	PC3
Applications	IHC
Host mice	Balb/C nu
Delivery	1 week
Guarantee	36 months from the date of manufacture
Background	These carefully prepared tissue slides offer researchers standardized and high-quality samples for antibody testing. Utilizing tumor models derived from immunocompromised mice, researchers can confidently assess antibody specificity and sensitivity, driving advancements in cancer diagnostics and therapeutics.
Yefei_Storage	Room temperature
Usage	Research use only

DIMABIO CONFIDENTIAL

