Cat. No. DME100192B



PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 bM192 Clone ID

Synonyme

Host Species Rabbit

Biotinylated Anti-ADAM9 antibody(DM192); Rabbit mAb Description 2-3 weeks

Delivery Uniprot ID Q13443 lgG type Rabbit IgG Clonality Monoclonal Reactivity Application ELISA; Flow Cyt

Recommended Dilutions ELISA 1:5000-10000; Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage & Shipping

specific instructions of reconstitution.

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing), tyophilized proteins are shipped at ambient temperature.

For example, the standard proteins are shipped at ambient temperature.

For example, the standard proteins are shipped at ambient temperature.

For example, the standard proteins are many analyses of this family are membrane-anchored proteins structurally related to stake venom disintegrins; and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions; including fertilization; muscle development; and neurogenesis. The protein encoded by this gene interacts with SH3 domain-containing proteins; binds mitotic arrest deficient 2 beta protein; and is also involved in TPA-induced ectodomain shedding of membrane-anchored heparin-binding EGF-like growth factor. Several alternatively spliced transcript variants have been identified for this gene.

Research use only

Email: info@dimabio.com Website: www.dimabio.com

Usage Coniugate

Background

, rotein s. All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer