

## **PRODUCT INFORMATION**

Warning: Undefined variable \$hasAttributeValueDescription in C:\www.root\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products-public.php on line 2806
\$MS-10.15MS 101.3MS 101 Common Name

Conjugate

B7H5;GI24;B7-H5;Dies1;PD-1H;SISP1;VISTA;PP2135;C10orf54;DD1alpha Synonyms

Applications Recommended Dilutions

DIMARIO CONTINUE DINAMENTO CONTINUE DE LA CONTINUE Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species

lgG1 IgG type Reactivity Target B7-H5 Uniprot ID

Biotinylated Anti-B7-H5(solnerstotug biosimilar) mAb Description

Delivery

6-8 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Usage

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

Email: info@dimabio.com Website: www.dimabio.com

