Cat. No. DME100092B



PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public:class-woocommerce-print-products-public-php on line 2806 bM92 Clone ID

RTN3A1

BTN3A1; BTF5; CD277; BTN3.1; BT3.1 Synonyme

Host Species Rabbit

Biotinylated Anti-BTN3A1 antibody(DM92); Rabbit mAb Description

Delivery 2-3 weeks Uniprot ID 000481 lgG type Rabbit IgG Clonality Monoclonal Reactivity Application ELISA; Flow Cyt

Recommended Dilutions ELISA 1:5000-10000; Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

The butyrophilin (BTN) genes are a group of major histocompatibility compiler (MHC)-associated genes that encode type I membrane proteins with 2 free that the control of the store of the st Background

Email: info@dimabio.com Website: www.dimabio.com

Usage Conjugate

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

