

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806 HIB210;HIB-210;HIB 210;MOR044254;MOR-044254;MOR 044254;MOR210;MOR-210;MOR 210;TJ 210;TJ-210;TJ210;TJ210001;WBP2191;WBP-2191;WBP 2191
Conjugate	Biotinylated
Synonyms	CSA;CSAR;CSR1;CD88
Applications	N/A
Recommended Dilutions	N/A
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	N/A
IgG type	IgG1
Reactivity	Human
Target	C5AR1
Uniprot ID	P21730
Description	Biotinylated Anti-C5AR1(izastobart biosimilar) mAb
Delivery	6-8 weeks
Yefei_Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

DIMABIO CONFIDENTIAL

