Cat. No. DMC100446B



PRODUCT INFORMATION

Warning: Undefined variable ShasAttributeValueDescription in C:\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public.class-woocommerce-print-products-public.php on line 2806 bMC446 Clone ID

CD117

C-Kit; CD117; MASTC; PBT; SCFR; KIT Synonyme

Host Species Rabbit

Biotinylated Anti-CD117 antibody(DMC446); IgG1 Chimeric mAb Description

Delivery 2-3 weeks Uniprot ID P10721

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Applications Flow Cyt Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

specific instructions of reconstitution.

Store at -20°C to -80°C for 12 months in jophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

The protein are shipped at ambient temperature.

The protein are shipped at ambient temperature, and the protein are shipped at ambient temperature.

The protein are shipped at ambient temperature, and the protein are shipped at ambient temperature.

The protein are shipped at ambient temperature are shipped at ambient temperature.

The protein are shipped at ambient temperature are shipped at ambient temperature, and the following temperature are shipped at ambient temperature.

The protein are shipped at a shipped at

Email: info@dimabio.com Website: www.dimabio.com

Usage Conjugate

Background

Research use only

Biotinylated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. **DIMA Disclaimer**



