Cat. No. DME100032B



PRODUCT INFORMATION

Warning: Undefined variable shasAttributeValueDescription in C\u00e4\u00fawwroot\u00fairor.dimabio.com\u00e4\u00fap-content\plugins\u00e4\u00faocommerce-print-products\u00e4\u00fap-billoc\u00e4\u00e Clone ID

CD123 IL3R; IL3RA; IL-3Ra; IL-3R-alpha; IL3RAY; IL3RX; IL3RY; CD123 antigen; CD123; hIL3Ra; hIL-3Ra; MGC34174; IL-3 R alpha Synonyme

Host Species

Biotinylated Anti-CD123 antibody(DM33); Rabbit mAb Description Delivery 2-3 weeks

P26951 Uniprot ID IgG type Rabbit IgG Clonality Monoclonal Reactivity Application ELISA; Flow Cyt

Recommended Dilutions ELISA 1:5000-10000; Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

Store at -20°C to .80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Interleukin 3 receptor alpha (low affinity) (IL3RA); also known as CD123 (Cluster of Differentiation 123) is a 70-kD glycoprotein member of the hematopoietin receptor superfamily. This protein associates with a beta subunit common to the receptors for IL-??

Research use only .xing or reverse

Background

Conjugate Biotinylated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

> Email: info@dimabio.com Website: www.dimabio.com

