Cat. No. DMC100217B



PRODUCT INFORMATION

Warning: Undefined variable shasAttributeValueDescription in Cs\www.root\wirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public.chass-woocommerce-print-products-public.php on line 2806 bMC217 Clone ID

CD24 CD24A Synonyme Host Species Rabbit

Biotinylated Anti-CD24 antibody(DMC217); IgG1 Chimeric mAb Description

Delivery 2-3 weeks Uniprot ID P25063

lgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Application ELISA; Flow Cyt

Recommended Dilutions ELISA 1:5000-10000; Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage & Shipping

Specific instructions or reconstitution.

Store at 2-90° to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at a maheint temperature.

This gene encodes a sidelity coprotein that is expressed on matter granulocytes and modulates growth and differentiation signals to these this gene expressed on second of the second of the

Email: info@dimabio.com Website: www.dimabio.com

Usage Conjugate Research use only Biotinylated

Background

Acten 5 All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

