

PRODUCT INFORMATION

|                              |                                                                                                                                                                                                                                                       |
|------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Common Name                  | <b>Warning:</b> Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimabio.com\wp-content\plugins\woocommerce-print-products\public\class-woocommerce-print-products-public.php on line 2806<br>SRF 617,SRF617                    |
| Conjugate                    | Biotinylated                                                                                                                                                                                                                                          |
| Synonyms                     | ENTPD1                                                                                                                                                                                                                                                |
| Applications                 | Flow Cyt                                                                                                                                                                                                                                              |
| Recommended Dilutions        | Flow Cyt 1:100                                                                                                                                                                                                                                        |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.                                                      |
| Host Species                 | Homo sapiens                                                                                                                                                                                                                                          |
| IgG type                     | IgG4                                                                                                                                                                                                                                                  |
| Reactivity                   | Human                                                                                                                                                                                                                                                 |
| Target                       | CD39                                                                                                                                                                                                                                                  |
| Uniprot ID                   | P49961                                                                                                                                                                                                                                                |
| Description                  | Biotinylated Anti-CD39(perenostobart biosimilar) mAb                                                                                                                                                                                                  |
| Delivery                     | 2-3 weeks                                                                                                                                                                                                                                             |
| Storage & Shipping           | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature. |
| Background                   | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.                                                                                                                                                 |
| Usage                        | Research use only                                                                                                                                                                                                                                     |

