Cat. No. DMC100685B



PRODUCT INFORMATION

Warning: Undefined variable shasAttributeValueDescription in Cs\www.root\wirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public.chass-woocommerce-print-products-public.php on line 2806
bMC685 Clone ID

CEACAM6 CD66c; CEAL; NCA Synonyme Rabbit

Host Species Biotinylated Anti-CEACAM6 antibody(DMC685); IgG1 Chimeric mAb Description

Delivery 2-3 weeks Uniprot ID P40199

lgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Applications Flow Cyt Recommended Dilutions Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage & Shipping

specific instructions of reconstitution.

Store at 20°C to 30°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing), Lyophilized proteins are shipped at ambient temperature repeated freezing and thawing), Lyophilized proteins are shipped at ambient temperature.

For example, the protein of the protei Background

> Email: info@dimabio.com Website: www.dimabio.com

Usage Coniugate

, rotein s All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

