

PRODUCT INFORMATION

Warning: Undefined variable \$hasAttributeValueDescription in Cs\www.root\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products\public.class-woocommerce-print-products-public.php on line 2806
bM179 Clone ID

CLDN18.2 Claudin 18.2 Synonyme Host Species Rabbit

Biotinylated Anti-CLDN18.2 antibody(DM179); Rabbit mAb Description

Delivery 2-3 weeks Uniprot ID P56856 IgG type Rabbit IgG Clonality Monoclonal Reactivity Human ELISA; Flow Cyt Applications

Recommended Dilutions ELISA 1:5000-10000; Flow Cyt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

> Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping

specific instructions or reconstitution.

Store at 2-0° to .80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature of the control of the

Usage Research use only

DIMA Disclaimer

biotinylated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infiningement.