

PRODUCT INFORMATION

Warning . Undefined variable \$hasAttributeValueDescription in Calwawroottmirror.dimablo.com/wp-content\plugins/woocommerce-print-products/public\class-woocommerce-print-products-public.php on line 2806 67E6 Clone ID

DLL3 SCD01 Synonyme Host Species Rabbit

Biotinylated Anti-DLL3 antibody(67E6), IgG1 Chimeric mAb Description

Delivery 2-3 weeks Uniprot ID Q9NYJ7

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Applications Flow Cyt Recommended Dilutions Flow Cyt 1/100

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

specinic instructions or reconstitutions. Store at 2-70°C to -80°C for 12°C to -80°C for 12°C for 12°C

Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only

Conjugate

Biotinylated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. Ain set. DIMA Disclaimer

