

PRODUCT INFORMATION

Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror.dimablo.com\wp-content\plugins\woocommerce-print-products-public.php on line 2806 LA12 Clone ID

Synonyme Host Species Rabbit

Biotinylated Anti-Dxd antibody(1A12); Rabbit mAb Description

Delivery 2-3 weeks Uniprot ID N.A. lgG type Rabbit IgG Clonality Monoclonal Reactivity ELISA Application Recommended Dilutions ELISA 1:5000-10000

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

specific instructions or recrostruction.

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature and the process of th

Usage Research use only Conjugate

Nescentro Use Unity
Biotinylated
All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

Email: info@dimabio.com Website: www.dimabio.com