

PRODUCT INFORMATION

Warning . Undefined variable \$hasAttributeValueDescription in Calwawroottmirror.dimablo.com/wp-content\plugins/woocommerce-print-products-public.php on line 2806 511 Clone ID

Feline IL31 nterleukin-31 IL-31 Synonyme

Host Species Rabbit

Biotinylated Anti-Feline IL31 antibody(5F11), Rabbit mAb Description

Delivery 2-3 weeks XP_011286140.1 Uniprot ID IgG type Rabbit IgG Clonality Monoclonal Reactivity Feline Application ELISA

Recommended Dilutions ELISA 1/5000-10000

Purification Purified from cell culture supernatant by affinity chromatography

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Storage & Shipping

specinic instructions or reconstitution. Store at 2-00° C. 0. 80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

1.31, which is made principally by activated Th2-type T cells, interacts with a heterodimeric receptor consisting of IL31RA (MIM 609510) and OSMR (MIM 601743) that is constitutively expressed on epithelial cells and keratinocytes. IL31 may be involved in the promotion of allergic skin disorders and in regulating other allergic diseases, such as a strain (Ollion et al., 2004 (PubMed 13184966)) supplied by OMIM, Mar 2008)

Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only Conjugate

Nescentro Use Unity
Biotinylated
All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

