

PRODUCT INFORMATION

Warning . Undefined variable \$hasAttributeValueDescription in Calwawroottmirror.dimablo.com/wp-content\plugins/woocommerce-print-products-public.php on line 2806 1945 Clone ID

FLT3

Flt-3;Flk-2;STK-1;CD135;FLK2;FLT-3 Synonyme

Host Species Rabbit

Biotinylated Anti-FLT3 antibody(19A5), IgG1 Chimeric mAb Description

Delivery 2-3 weeks Uniprot ID P36888

lgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Application Flow Cyt Recommended Dilutions Flow Cyt 1/100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage & Shipping

Email: info@dimabio.com Website: www.dimabio.com

specific instructions of reconstitution.

Store at 2.0º Cto 3.0º Cto 1.0º Cto 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -8.0°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

FIL-3(Receptor-type tyrosine-protein kinase FLT3) is also known as FLK-2(Fetal liver kinase-2), STK-1(Stem cell tyrosine kinase 1), CD135. FLT3 is a FLT3 is a classification within the protein kinase being to the receptor tyrosine kinase that set as cell-surface receptor for the cytokine receptor which belongs to the receptor tyrosine kinase class III. Tyrosine-protein kinase that acts as cell-surface receptor for the cytokine FLT3LG and regulates differentiation, proliferation and survival of hematopoietic progenitor cells and of dendritic cells. Promotes schosphorylation of Edwistream Kinases SYCLC and Art1, and activation of the Counstream effector MT0/R. Promotes activation of RAS signaling and phosphorylation of downstream kinases of the standard signaling signaling and prophers of the standard signaling pathways.

Research use only

Usage

Biotinylated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

