

PRODUCT INFORMATION

Warning . Undefined variable \$hasAttributeValueDescription in Calwawroottmirror.dimablo.com/wp-content\plugins/woocommerce-print-products-public.php on line 2806
1843 Clone ID

FOLR1 FBP;FOLR;FRalpha Synonyme

Host Species Rabbit

Biotinylated Anti-FOLR1 antibody(18A3), IgG1 Chimeric mAb Description

Delivery 2-3 weeks Uniprot ID P15328

IgG type Rabbit/Human Fc chimeric IgG1

Clonality Monoclonal Reactivity Human Application WB Recommended Dilutions WB 1:1000

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage & Shipping

specific instructions of reconstitution.

Store at 20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature temperature repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature temperature and the store of the stor Background

Usage Coniugate

/ protein seq All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. DIMA Disclaimer

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

Email: info@dimabio.com Website: www.dimabio.com

