

PRODUCT INFORMATION

Common Name	Warning: Undefined variable \$hasAttributeValueDescription in C:\wwwroot\mirror\dimabio.com\wp-content\plugins\woocommerce-print-products\publicclass-woocommerce-print-products-public.php on line 2806 PF-06946860
Conjugate	Biotinylated
Synonyms	GDF-15;MIC-1;MIC1;NAG-1;PDF;PLAB;PTGFB
Applications	ELISA
Recommended Dilutions	ELISA 1:5000-10000
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Humanized
IgG type	IgG1
Reactivity	Human
Target	GDF15
Uniprot ID	Q99988
Description	Biotinylated Anti-GDF15(ponsegromab biosimilar) mAb
Delivery	In Stock
Yefei Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only

Biotinylated Anti-GDF15(ponsegromab biosimilar) mAb ELISA

0.2 µg of Human GDF15, hFc tagged protein per well




Figure 1. ELISA plate pre-coated by 2 µg/mL (100 µL/well) Human GDF15 Protein, hFc Tag (PME100628) can bind Biotinylated Anti-GDF15(ponsegromab biosimilar) mAb (BME100090B) in a linear range of 8.0-40 ng/mL.

